


KatKath™

ADJUSTABLE URINARY CATHETER

Be kind to your patients



Features

- Non-reactive material
- Suture groove on hub
- Low coefficient of friction
- Kink resistant
- Crush resistant
- Lightweight
- Radiopaque
- No stylet required

Benefits

- Easy placement
- Adjustable length
- Flexible for maximum comfort
- Thin walled - larger bore lumen
- Great longitudinal strength
- Rounded atraumatic tip

Available in boxes of 10 units

| Code | Description | Box quantity |
|-------------|--|--------------|
| VET98AK3518 | Adjustable urinary catheter 3.5Fr/18cm | 10 |
| VET98KK3511 | Tom cat urinary catheter 3.5Fr/11cm | 10 |
| VET98KK3514 | Tom cat urinary catheter 3.5Fr/14cm | 10 |

KatKath™

ADJUSTABLE URINARY CATHETER



Be kind to your patients

FULLY ADJUSTABLE SUTURE WING

The Adjustable KatKath™ catheter allows adjustment through the full range of lengths as needed for each individual patient from the very smallest to the largest of cats. Sutures should be tied around the two suture grooves on the body of the suture wing component thus securing the catheter at the required length.

MORE FLEXIBLE

The catheter is made from a unique non-reactive material which has similar properties to Teflon® but is more flexible. More flexibility allows greater comfort.

GREAT LONGITUDINAL STRENGTH

The KatKath's increased flexibility does not diminish its ability in the effective treatment of blocked cats due to its great longitudinal strength and will not concertina even when encountering gross resistance.

CRUSH RESISTANT

The KatKath™ is almost impossible to crush and is also highly resistant to kinking. This allows placement without the use of a stylet. (Useful if flushing is required during the act of placement).

WELL TOLERATED

The KatKath™ may be left in place for hours or days, dependant on individual needs, and is extremely well tolerated by even the most difficult patients.

KATKATH™ SAVES TIME AND MONEY

There is now no longer any need to purchase two separate catheters for flushing blocked cats and leaving a different catheter in place. The KatKath™ is highly effective in both situations saving time, money and effort, by cutting two procedures down to one.

To place your order contact us at: 01793 748900